



Neuroendoskopie  
Neuroendoscopy  
Neuroendoscopie  
Neuroendoscopia  
Neuroendoscopia

## Vertebroskop

Vertebroscope  
Vertebroscope  
Vertebroscopio  
Vertebroscopio

Endoskope Ø 4,0 mm  
Endoscope Ø 4.0 mm  
Endoscope Ø 4.0 mm  
Endoscopio Ø 4.0 mm  
Endoscopio Ø 4.0 mm



### 91.876.30

30°	Blickrichtung / Direction of view / Direction de vue / Dirección de vista / Direzione di vista
70 mm	Arbeitslänge / Working length / Longueur de travail / Longitud de Trabajo / Lunghezza di lavoro
90°	Bildwinkel / Field of view / Champ de vision / Campo de visión / Campo di vista



### 91.876.31

Arbeitskanal und Haltehülse für **91.876.30**, erlaubt einen Drehbereich von 360°  
Sheath and Tube for **91.876.30**, rotatable 360°  
Tenue de disque rotatif, de tenir l'endoscope **91.876.30**, tournante 360°  
Soporte giratorio para optica **91.876.30**, giratorio 360°  
Sopporto girevole per ottica **91.876.30**, girevole 360°